

SBCNY

NIGMS funded Center

Experimental Methods in Systems Biology

Part of the Coursera Certificate in Systems Biology

Marc Birtwistle, PhD

Department of Pharmacology & Systems Therapeutics

Fall 2014, Week 3 Mass Spectrometry-Based
Proteomics



Icahn School
of Medicine at
**Mount
Sinai**




Quantification

MS-based protein quantification

- Stable isotope-labeling methods
 - SILAC: Stable isotope labeling with amino acids in cell culture
 - iTRAQ: Isobaric Tags for Relative and Absolute Quantification
- Selective Reaction Monitoring
 - Need to know what you're looking for and typically employs internal standards
 - We're focusing more on omic-level proteomics so we won't cover this
- Label-free method
 - Spectral counting—counts of peptides that align to a protein, analogous to aligned reads in mRNA seq
 - iBAQ—intensity based absolute quantification


iTRAQ reagents

(Isobaric Tags for Relative and Absolute Quantification)


	Reporter Group	Balance Group	Amine-reactive Group
			
1	^{13}C	114.1	$^{13}\text{C}, ^{18}\text{O}$ 31.0
2	$^{13}\text{C}_2$	115.1	^{18}O 30.0
3	$^{13}\text{C}_2, ^{15}\text{N}$	116.1	^{13}C 29.0
4	$^{13}\text{C}_3, ^{15}\text{N}$	117.1	28.0

- Allows multiplexing
 - 4 or 8
- Labeled peptides co-elute on LC separation
- Differences in mass only apparent after fragmentation

SILAC—Stable Isotope Labeling with Amino Acids in Cell Culture




Geiger et al., Nat Protocols, 2011



- Heavy isotopes of lysines and arginine residues are incorporated into proteome
- The mass shift allows relative quantification


Spike-in SILAC



Geiger et al., Nat Protocols, 2011

- Some samples very expensive and/or not possible to label
 - E.g. mice, human tissues
- The spike-in method of labeled standards can solve that
- Not all tissues have similar expression profiles as a single standard
 - There is a “super SILAC standard” that addresses this issue by mixing many standards (Geiger et al., 2010, Nat Methods)

iBAQ—Intensity-Based Absolute Quantification



Schwanhausser et al., 2011, Nature

- Sum all identified peptide intensities
 - Maximum detector peak intensities of the peptide elution profile
- Divide by the number of theoretically observable peptides
 - Similar to mRNA seq divide by transcript length
- Seems to perform (slightly) better than other label-free methods (Arike et al., J Proteomics, 2012)
- Can be converted into absolute copy number / cell by keeping track of how much protein is kept from lysate until mass spec or with known spike-ins