

SBCNY

NIGMS funded Center

Experimental Methods in Systems Biology

Part of the Coursera Certificate in Systems Biology

Marc Birtwistle, PhD

Department of Pharmacology & Systems Therapeutics

Fall 2014, Week 6, Live-Cell Imaging



Icahn School
of Medicine at
**Mount
Sinai**


Outline

- Fluorescence Microscopy—Instrumentation
 - Components of a fluorescence microscope
 - Keeping cells “happy”—environmental control
 - Types of imaging
 - Acquisition properties
- Fluorescent Probes
 - Small molecules
 - Genetically-encoded sensors
- Quantification
- A note about
 - super resolution—relatively new technique that allows imaging past the diffraction limit but we won’t cover it here.
 - single molecule imaging—it is possible with fluorescence microscopy in live cells but we will not go into detail on it. Useful for modeling stochastic phenomena.
 - Fluorescence recovery after photobleaching (FRAP)—usually done in live cells to estimate diffusivity of tagged proteins but typically not repeatedly...we won’t cover it here

Fluorescence Microscopy— Instrumentation

General Components of a Fluorescence Microscope

http://en.wikipedia.org/wiki/Fluorescence_microscope#mediaviewer/File:FluorescenceFilters_2008-09-28.svg




- Light source
 - Lasers
 - Defined wavelengths
 - “White” lasers
 - Broad coverage lamps
 - Xenon
 - Mercury
 - LED
 - Defined wavelengths
 - “White” LEDs
- Excitation filter
 - Allows only certain wavelengths of light through
 - Bandpass or longpass
- Dichroic Mirror
 - reflects excitation light but lets through emission light
 - Can also be polychroic
- Emission filter
 - Similar to excitation filter, but lets a different set of wavelengths through
- The excitation filter, dichroic mirror and emission filter are often contained together in a filter cube. They can be separate though.
- Objective
 - Lens that magnifies the image of the specimen
- Stage
 - Where the specimen sits, can be fixed or movable by hand or motorized
- Ocular
 - The eyepiece followed by the detector being one’s eye
 - Can also be a camera
- This schematic is an upright microscope, but for live cell imaging inverted microscopes are typically used (objective below specimen)

- Successful live-cell imaging experiments depend on the ability to maintain the cells in a healthy state and functioning normally on the microscope stage while being illuminated in the presence of synthetic fluorophores and/or fluorescent proteins.


Environment: Keeping cells alive

- Media (pH, osmolarity)
- Temperature
- Atmospheric conditions (gas mixture and humidity)

Incubation chambers can maintain cells under desired conditions for long periods of time.



Microscope Configuration for Live-Cell Imaging



<http://zeiss-campus.magnet.fsu.edu/articles/livecellimaging/imagingsystems.html>

Vibration Isolation Systems

Figure 4




(a)

Neoprene Damping Pads



(b)

Sorbothane Damping Pads



(c)

Vibration Isolation Table

<http://www.microscopyu.com/articles/livecellimaging/livecellmaintenance.html>