

SBCNY

NIGMS funded Center

Experimental Methods in Systems Biology

Part of the Coursera Certificate in Systems Biology

Marc Birtwistle, PhD

Department of Pharmacology & Systems Therapeutics

Fall 2014, Week 5 Flow and Mass Cytometry




Icahn School
of Medicine at
**Mount
Sinai**

Outline

- Flow Cytometry
 - How does it work?
 - What can you measure?
- Mass Cytometry
 - How does it work?
 - What can you measure?
- Analysis
 - Gating strategies
 - Debris vs. cells
 - Singlets vs. doublets
 - Automated gating strategies
 - Compensation—what to do when signals overlap
 - Quantification
 - DNA stains
 - Antibody stains

Flow Cytometry



- Flow cytometry is a technology for characterizing and sorting cells as they move in a fluid stream through laser beams
- Allows simultaneous quantitative measurements of multiple parameters on a single cell level in a heterogeneous cell population
- Allows quantitative measurements of cell frequencies

Flow Cytometers: Analyzers

BD FACS Calibur
2 lasers; 4 colors



BD FACS Canto
2 lasers; 6 colors
3 lasers; 8 colors



BD LSRII
5 lasers;
18 colors



BD LSR
Fortessa
5 lasers;
18 colors



Mass Cytometer

DVS CyTOF2
No lasers; > 40 colors



Flow Cytometers: Cell sorters




BD FACSAria II

3L, 4L & 5L

Up to 4-way sorting

Non-BSL and BSL2+



BD Influx


5L

Up to 6-way sorting

Non-BSL

Flow Cytometer


Components & Principles




Liquid management: fluidics cart



The Flow Cell: hydrodynamic focusing and interrogation




Fluorescence Activated Cell Sorting (FACS)




Light Scatter: FSC & SSC


Side Scatter Histogram







Forward Scatter





Fluorescence Detection: Filters



Copyright ©2011 Life Technologies Corporation. Used under permission.



Example: Fluorescein (FITC)



Linear vs. Logarithmic scales

Generally speaking...



- FSC and SSC parameters are displayed on a linear scale
- Fluorescent parameters are displayed on a log scale

Some notable exceptions:

- RBC/platelet analysis is often done using FSC on log scale
- DNA content should be examined on a linear scale


Linear: $x = y$

Log: $y = \ln(x)$




What Can Be Measured with a Flow Cytometer

Events




2D scatter plot of FSC vs SSC




Viability and Apoptosis

- Many fluorescent dyes are cell impermeant unless the membrane is compromised
 - E.g. propidium iodide, Draq7
- This allows detection of viability
- That can be combined with fluorophore-labeled Annexin-V
 - binds to exposed phosphatidylserine, a marker of apoptosis




Cell Cycle Analysis



- There are many fluorescent dyes that can stoichiometrically bind to DNA
 - E.g. propidium iodide
- DNA dyes can be combined with BrdU staining to separate cells actively replicating their DNA
 - A pulse of BrdU gets incorporated into the DNA, which can then be recognized by an antibody

Immunophenotyping



http://docs.abcam.com/pdf/immunology/lymphoid_cell.pdf

- One can stain cells with labeled antibodies against lineage-specific extracellular markers


Intracellular Immunostaining

SCIENCE VOL 321 22 AUGUST 2008

Variability and Robustness in T Cell Activation from Regulated Heterogeneity in Protein Levels

Ofer Feinerman,^{1*} Joël Veiga,^{1*} Jeffrey R. Dorfman,^{1†}
Ronald N. Germain,² Grégoire Altan-Bonnet^{1‡}

- After cell fixation, one can detect intracellular epitopes
- E.g. panels of phosphorylated proteins



Fluorescent Cell Barcoding (FCB) for Multiplex Flow Cytometry

